


AVEX

SURGICAL

Shaped with precision



EXTRACTION FORCEPS (ENGLISH PATTERN)



Zahnärztliche Instruments

Dentales Extracción Instrumentos

Instruments de Dentaire Extraction

Dentale Estrazione Strumenti

APPLICATION CHART FOR EXTRACTION FORCEPS





Fig	Application
1	Upper centrals, canines and roots
2	Upper laterals
7	Upper premolars and roots
13	Lower premolars
17	Upper molars (Right)
18	Upper molars (Left)
19	Upper wisdoms
22	Lower molars (Hawk's bill)
29	Upper roots
29S	Upper roots
30	Upper roots
37	Children's upper incisors and canines
39	Children's upper molars
44N	Upper roots
51	Upper wisdoms (Bayonet)
67	Upper wisdoms
73	Lower molars (Hawk's bill)
73S	Lower molars (Hawk's bill)
74	Lower roots
74N	Small lower roots and anteriors
75	Lower premolars
76	Upper roots (Read)
76S	Upper roots (Read)
76N	Small upper roots (Read)

Fig	Application
79	Lower wisdoms
86	Lower molars (Cowhorn)
87	Lower molars (Cowhorn)
89	Upper molars (Right)(Upper Cowhorns)
90	Upper molars (Left)(Upper Cowhorns)
94	Upper molars (Right)
95	Upper molars (Left)
101	Upper molars and roots
107	Upper canines
113	Upper roots (Martin)
123	Children's lower anteriors, premolars, roots
136	Upper premolars and roots (Guy's)
137	Upper anteriors and roots (Guy's)
138	Children's upper anteriors, premolars, roots
147	Upper roots (Gibb)
157	Children's upper molars
158	Children's upper molars
159	Children's upper premolars
160	Children's lower molars (Hawk's bill)
161	Children's lower molars (Hawk's bill)
162	Children's lower anteriors and roots
163	Children's upper anteriors
3	Rongeurs forceps

Shaped with precision

AVEX
SURGICAL

HANDLE OPTIONS

SMALL MODEL



LARGE MODEL



FINGER-SHAPED



DIAMOND TIPPED






450-01 Fig. 01
Upper centrals and canines, wide



450-02 Fig. 02
Upper centrals and canines,



450-02 Fig. 03
Upper centrals and canines, curved



450-04 Fig. 04
Lower incisors and canines



450-07 Fig. 07
Upper Premolars



450-08 Fig. 08
Lower premolars either side



450-13 Fig. 13
Lower premolars either side



450-13S Fig. 13S
Lower premolars, for children




450-17 Fig. 17
Upper molars, right



450-18 **Fig. 18**
Upper molars left



450-18A **Fig. 18A**
Upper molars either side



450-18B **Fig. 18B**
Upper molars



450-19 **Fig. 19**
Upper third molars & upper
wisdoms either side



450-20 **Fig. 20**
Lower third molars



450-21 **Fig. 21**
Lower molars



450-22 **Fig. 22**
Lower molars



450-22S **Fig. 22S**
Lower molars, for children



450-23 **Fig. 23**
Lower molars, right



450-24 **Fig. 24**
Lower molars, left



450-29 **Fig. 29**
Upper roots and incisors



450-29N **Fig. 29N**
Root of Upper Incisors



450-29S **Fig. 29S**
Upper roots, for children



450-30 **Fig. 30**
Upper roots



450-31 **Fig. 31**
Lower roots & incisors with
Serrated tips




450-32 **Fig. 32**
Lower molars



450-33 **Fig. 33**
Lower roots



450-33A **Fig. 33A**
Lower roots



450-33C **Fig. 33C**
Lower roots




450-33L **Fig. 33L**
Lower roots



450-33M **Fig. 33M**
Lower roots with long beak



450-34 **Fig. 34**
Incisors and Bicuspid, deep gripping



450-34N **Fig. 34N**
Incisors and Bicuspid, deep gripping



450-35 **Fig. 35**
Incisors and Bicuspid, deep gripping



450-35N **Fig. 35N**
Incisors and Bicuspid, deep gripping



450-36 **Fig. 36**
Incisors and Bicuspid, deep gripping



450-37 **Fig. 37**
Upper incisors, for children




450-38 Fig. 38
Lower incisors and
canines, for children




450-39 Fig. 39
Upper molars, for children




450-39A Fig. 39A
Upper molars, for children



450-39L Fig. 39L
Upper molars left, for children



450-39R Fig. 39R
Upper molars right, for children



450-40 Fig. 40
Lower molars, for children



450-41 Fig. 41
Upper roots



450-44 Fig. 44
Upper roots



450-44F Fig. 44F
Upper roots




450-45 **Fig. 45**
Lower roots



450-45F **Fig. 45F**
Lower roots



450-46 **Fig. 46**
Lower roots



450-46L **Fig. 46L**
Lower roots



450-46LX **Fig. 46LX**
for very fine roots



450-46N **Fig. 46N**
Bicuspid, deep gripping




450-49 **Fig. 49**
Upper roots



450-50 **Fig. 50**
Upper roots




450-51 **Fig. 51**
Upper roots




450-51A **Fig. 51A**
Upper roots




450-51C **Fig. 51C**
Upper roots



450-51L **Fig. 51L**
Upper roots



450-51M **Fig. 51M**
Especially long beaks



450-52 **Fig. 52**
Upper roots



450-52A **Fig. 52A**
Upper teeth & Resection




450-53 **Fig. 53**
Upper roots




450-54 **Fig. 54**
For separating upper molars




450-55 **Fig. 55**
For separating upper molars



450-56 **Fig. 56**
For separating lower molars



450-65L **Fig. 65L**
Upper molars, left



450-65R **Fig. 65R**
Upper molars, right




450-67 **Fig. 67**
Upper third molars and
upper wisdoms either side



450-67A **Fig. 67A**
Upper third molars and upper
wisdoms either side



450-67LX **Fig. 67LX**
Upper third molar, left



450-67RX **Fig. 67RX**
Upper third molar, right



450-67N **Fig. 67N**
Third molar, deep-gripping



450-68 **Fig. 68**
Lower roots

Shaped with precision

AVEX

SURGICAL



450-73 **Fig. 73**
Lower molars



450-73L **Fig. 73L**
Lower molars



450-73S **Fig. 73S**
Lower molars



450-74 **Fig. 74**
Lower roots and incisors




450-74D **Fig. 74D**
Lower roots



450-74M **Fig. 74M**
Lower roots and incisors



450-74N **Fig. 74N**
Lower roots and incisors



450-75 **Fig. 75**
Lower incisors and premolars



450-75A **Fig. 75A**
Lower incisors and premolars




450-76 **Fig. 76**
Upper roots




450-76N **Fig. 76N**
Upper roots



450-76S **Fig. 76S**
Upper roots



450-77 **Fig. 77**
Upper roots



450-79 **Fig. 79**
Lower third molars and
lower wisdoms either side



450-79A **Fig. 79A**
Lower third molars



450-79N **Fig. 79N**
Lower Wisdoms



450-86 **Fig. 86**
Lower molars decayed
or broken down crowns



450-86B **Fig. 86B**
Lower molars



450-86C Fig. 86C
Lower molars decayed or
broken down crowns



450-87 Fig. 87
Lower molars decayed or
broken down crowns




450-89 Fig. 89
Upper molars, right



450-90 Fig. 90
Upper molars, left




450-91 Fig. 91
For separating lower molars



450-92 Fig. 92
Upper molars, right



450-93 Fig. 93
Upper molars, left



450-94 Fig. 94
Upper molars, right




450-95 Fig. 95
Upper molars, left



450-96 **Fig. 96**
Upper molars, right




450-97 **Fig. 97**
Upper molars, left




450-99 1/4 **Fig. 99 1/4**
Lower molars



450-99 **Fig. 99**
Lower molars




450-101 **Fig. 101**
Upper roots




450-101A **Fig. 101A**
Upper roots




450-107 **Fig. 107**
Upper Canines




450-108 **Fig. 108**
Upper roots



450-110 **Fig. 110**
Upper roots



450-111 **Fig. 111**
Upper roots



450-112 **Fig. 112**
Upper roots



450-113 **Fig. 113**
Upper roots



450-114 **Fig. 114**
Lower roots & centrals



450-121 **Fig. 121**
Upper wisdoms



450-122 **Fig. 122**
Lower centrals & roots



450-123 **Fig. 123**
Lower teeth & roots,
for children




450-130 **Fig. 130**
Upper wisdoms 190 mm



450-133 **Fig. 133**
for very fine roots



450-136 **Fig. 136**
Upper Centrals,
Premolars & roots



450-136H **Fig. 136H**
Upper roots



450-137 **Fig. 137**
Lower centrals & roots



450-138 **Fig. 138**
Upper teeth & roots,
for children



450-139 **Fig. 139**
Upper roots



450-143 **Fig. 143**
Lower roots & centrals



450-144 **Fig. 144**
Lower roots & centrals



450-145 **Fig. 145**
Lower roots




450-146 **Fig. 146**
Root fragments,
deep gripping




450-147 **Fig. 147**
Upper roots




450-148 **Fig. 148**
Lower molars



450-151 **Fig. 151**
for very fine roots



450-153 **Fig. 153**
Root fragments,
deep gripping



450-155 **Fig. 155**
Lower roots



450-157 **Fig. 157**
Upper molars,
for children



450-158 **Fig. 158**
Upper premolars & molars,
for children



450-159 **Fig. 159**
Upper premolars,
for children



450-160 **Fig. 160**
Lower molars,
for children



450-161 **Fig. 161**
lower molars, for children



450-162 **Fig. 162**
Lower roots & centrals,
for children



450-163 **Fig. 163**
Upper centrals, for children



450-166 **Fig. 166**
Upper universal application



450-167 **Fig. 167**
Lower teeth, universal pattern



450-168 **Fig. 168**
Upper central & lateral teeth




450-169 **Fig. 169**
Premolars, lower canines
& all roots



450-171 **Fig. 171**
Wisdoms & molar roots



450-173 **Fig. 173**
Lingually inclined roots and
wisdom teeth



450-180 Fig. 180
Special for upper centrals and
canines deep gripping



450-181 Fig. 181
Special for upper premolars,
deep gripping



450-182 Fig. 182
Lower centrals canines and
premolars deep gripping



450-185 Fig. 185
Lower roots



450-187 Fig. 187
Lower roots fragments



450-189 Fig. 189
Fragments upper roots



450-190 Fig. 190
Fragments upper roots



450-221 Fig. 221
Universal upper roots 150 mm



450-223 Fig. 223
Universal for low roots



450-350 **Fig. 350**
Lower molars & wisdoms left



450-351 **Fig. 351**
Lower molars & wisdoms right

AVEX SURGICAL

AVEX SURGICAL INTERNATIONAL

DEFENCE ROAD, SIALKOT 51310, PAKISTAN

TEL: +92-328-1584072 / +92-305-4813441

Email: ask@avexsurgical.com / info@avexsurgical.com

Website: <http://www.avexsurgical.com>



Certificates & Memberships



Surgical Instruments • Chirurgische Instrumente • Instrumentos para quirúrgicos

Instruments de chirurgie • Strumenti chirurgici